


## (Z)-2,3-Dehydroadipoyl-CoA

## Chemical Properties

CAS No. :	138149-18-5
Formula:	C <sub>27</sub> H <sub>42</sub> N <sub>7</sub> O <sub>19</sub> P <sub>3</sub> S
Molecular Weight:	893.64
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(Z)-2,3-Dehydroadipoyl-CoA ((Z)-2,3-Dehydroadipoyl-coenzyme A) is a type of cis-2-enoyl coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.119 mL	5.5951 mL	11.1902 mL
5 mM	0.2238 mL	1.119 mL	2.238 mL
10 mM	0.1119 mL	0.5595 mL	1.119 mL
50 mM	0.0224 mL	0.1119 mL	0.2238 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481