

Flt-3L-Ig (hum/hum)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Flt-3L-Ig (hum/hum) (hFlt3L) is a Flt3 activator. It enhances the release of inflammatory cytokines from myeloid and dendritic cells in BRGSF-CBC mice induced by OKT3, as well as affects the quantity and distribution of human immune cells, exhibiting cytokine release syndrome (CRS).
Targets(IC50)	FLT

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481