

2-trans-9,12,15,18-all-cis-Tetracosapentaenoyl-CoA

Chemical Properties

CAS No. :

Formula: C45H72N7O17P3S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-trans-9,12,15,18-all-cis-Tetracosapentaenoyl-CoA (2-trans-9,12,15,18-all-cis-Tetracosapentaenoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481