


Phenylmethyl-D7 bromide

Chemical Properties

CAS No. : 35656-93-0
 Formula: C7D7Br
 Molecular Weight: 178.08
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenylmethyl-D7 bromide is the deuterium-labeled form of phenylmethyl bromide.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6155 mL	28.0773 mL	56.1545 mL
5 mM	1.1231 mL	5.6155 mL	11.2309 mL
10 mM	0.5615 mL	2.8077 mL	5.6155 mL
50 mM	0.1123 mL	0.5615 mL	1.1231 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481