

## DDR1 ligand 1

### Chemical Properties

CAS No. :

Formula: C<sub>22</sub>H<sub>16</sub>F<sub>3</sub>N<sub>5</sub>O<sub>3</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	DDR1ligand 1 is a ligand used as a target protein for PROTACs (Ligands for Target Protein for PROTACs). It is utilized in the synthesis of PROTAC DDR1 degrader-1.
Targets(IC50)	Ligands for Target Protein for PROTAC

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481