

## Cy3 Dextran (MW 20000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cy3Dextran (MW 20000) is a fluorescent labeling reagent that combines the Cy3 fluorescent dye with Dextran. The Cy3 fluorophore is commonly utilized in applications such as immunolabeling, nucleic acid labeling, fluorescence microscopy, and flow cytometry. Dextran inhibits platelet aggregation and clotting factors, serving as a plasma volume expander (Ex/Em = 550/570 nm).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481