

STAT3-IN-42

## Chemical Properties

CAS No. :

Formula: C34H42N7O11PS

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	STAT3-IN-42 is a ligand for target proteins in PROTACs (Proteolysis Targeting Chimeras). It is utilized in the synthesis of SD-2301.
Targets(IC50)	Ligands for Target Protein for PROTAC

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481