

Poly(methyl vinyl ether-alt-maleic anhydride) (MW 1080000)

Chemical Properties

CAS No. :

Formula: (C₄H₂O₃.C₃H₆O)_n

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Poly(methyl vinyl ether-alt-maleic anhydride) (MW 1080000) is suitable for preparing antibacterial hydrogels that exhibit excellent mechanical strength and swelling capacity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481