

## E3 Ligase Ligand-linker Conjugate 184

## Chemical Properties

CAS No. :

Formula: C<sub>29</sub>H<sub>41</sub>N<sub>5</sub>O<sub>5</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	E3LigaseLigand-linker Conjugate 184 is a conjugate consisting of an E3 ligase ligand and a linker (E3LigaseLigand-Linker Conjugates), composed of CDK2ligand-2 and its corresponding linker. E3LigaseLigand-linker Conjugate 184 can act as a Cereblon ligand to recruit CRBN protein and serves as a key intermediate in the synthesis of CDK2 degrader2.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481