


Parvodicin C1

Chemical Properties

CAS No. :	110882-84-3
Formula:	C83H88Cl2N8O29
Molecular Weight:	1732.53
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Parvodicin C1 is an intermediate of Dalbavancin.
Targets(IC50)	Antibacterial, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5772 mL	2.886 mL	5.7719 mL
5 mM	0.1154 mL	0.5772 mL	1.1544 mL
10 mM	0.0577 mL	0.2886 mL	0.5772 mL
50 mM	0.0115 mL	0.0577 mL	0.1154 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481