

ADP- β -S trisodium

Chemical Properties

CAS No. :

Formula: C₁₀H₁₂N₅Na₃O₉P₂S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ADP- β -S (trisodium) is an analog of ADP and acts as a partial agonist. It can induce human platelet aggregation and inhibit PGE ₁ -stimulated adenylate cyclase.
Targets(IC ₅₀)	PGE Synthase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481