

## Heparin-Cy7

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

|               |   |
|---------------|---|
| Description   | Heparin-Cy7 is a fluorescent dye composed of heparin and CY7. It is utilized to investigate the localization, distribution, and interactions of heparin with other biomolecules within biological systems (Ex/Em = 740/770 nm). |
| Targets(IC50) | Others  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481