

TLCD1-IN-1

## Chemical Properties

CAS No. :

Formula: C<sub>20</sub>H<sub>19</sub>ClNO<sub>5</sub>PS<sub>2</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | TLCD1-IN-1 (compound MB-4A) is a covalent inhibitor of TLCD1. It demonstrates stereoselective binding and activity inhibition of TLCD1 in HepG2 cells. |
| Targets(IC50) | Endogenous Metabolite  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481