


## 3-Decenoyl-CoA

## Chemical Properties

CAS No. : 6410-54-4  
 Formula: C<sub>31</sub>H<sub>52</sub>N<sub>7</sub>O<sub>17</sub>P<sub>3</sub>S  
 Molecular Weight: 919.77  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

|               |   |
|---------------|---|
| Description   | 3-Decenoyl-CoA (3-Decenoyl-coenzyme A) is a derivative of coenzyme A. |
| Targets(IC50) | Others  |

## Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg       |
|-------|-----------|-----------|------------|
| 1 mM  | 1.0872 mL | 5.4361 mL | 10.8723 mL |
| 5 mM  | 0.2174 mL | 1.0872 mL | 2.1745 mL  |
| 10 mM | 0.1087 mL | 0.5436 mL | 1.0872 mL  |
| 50 mM | 0.0217 mL | 0.1087 mL | 0.2174 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481