

08:0 PI ammonium

Chemical Properties


CAS No. : 2260795-76-2

Formula: C₂₅H₅₀NO₁₃P

Molecular Weight: 603.64

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	08:0 PI ammonium is a glycerophospholipid featuring a head group that carries a cyclohexanol moiety.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6566 mL	8.2831 mL	16.5662 mL
5 mM	0.3313 mL	1.6566 mL	3.3132 mL
10 mM	0.1657 mL	0.8283 mL	1.6566 mL
50 mM	0.0331 mL	0.1657 mL	0.3313 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481