

E3 Ligase Ligand-linker Conjugate 177

Chemical Properties

CAS No. :

Formula: C₂₂H₂₈ClN₅O₄

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3LigaseLigand-linker Conjugate 177 is a compound that functions as an E3 ligase ligand-linker conjugate. It is utilized in the synthesis of the STING-PROTAC degrader, ST9.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481