


SD-2301

Chemical Properties

CAS No. : 3054116-24-1
Formula: C₆₄H₈₀N₁₃O₁₅PS₂
Molecular Weight: 1366.50
Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	SD-2301 is a PROTAC degrader of STAT3.
Targets(IC50)	STAT,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7318 mL	3.659 mL	7.318 mL
5 mM	0.1464 mL	0.7318 mL	1.4636 mL
10 mM	0.0732 mL	0.3659 mL	0.7318 mL
50 mM	0.0146 mL	0.0732 mL	0.1464 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481