

FITC-PEG-CL (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	FITC-PEG-CL (MW 10000) is a fluorescent dye composed of FITC, PEG, and chlorine (-Cl). The chlorine group acts as an electrophilic moiety suitable for nucleophilic substitution reactions, facilitating subsequent chemical modifications (Ex/Em = 488/525 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481