

Cy5.5 HA MA (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5.5HA MA (MW 10000) is a methacrylate derivative of hyaluronic acid labeled with CY5.5, and is extensively utilized in fields such as biomarking and imaging, drug delivery, and tissue engineering and regenerative medicine (Ex/Em = 680/710 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481