


2-Carboxytetracosanoyl-CoA

Chemical Properties

CAS No. : 57514-08-6
Formula: C₄₆H₈₂N₇O₁₉P₃S
Molecular Weight: 1162.17
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Carboxytetracosanoyl-CoA (2-Carboxytetracosanoyl-coenzyme A) is a type of 2-carboxyacyl coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8605 mL	4.3023 mL	8.6046 mL
5 mM	0.1721 mL	0.8605 mL	1.7209 mL
10 mM	0.086 mL	0.4302 mL	0.8605 mL
50 mM	0.0172 mL	0.086 mL	0.1721 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481