

## KRAS inhibitor-42

### Chemical Properties

CAS No. :

Formula: C<sub>34</sub>H<sub>47</sub>ClN<sub>8</sub>O<sub>4</sub>S<sub>2</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	KRAS inhibitor-42 (compound 8) is an effective USP7 inhibitor that exhibits high affinity for GDP-bound KRASG12D, with a $K_i$ value of 2.7 $\mu$ M.
Targets(IC50)	Ras

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481