

3-(Benzoylamino)-L-Alanine-CoA

Chemical Properties

CAS No. :

Formula: C₃₁H₄₆N₉O₁₈P₃S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | 3-(Benzoylamino)-L-Alanine-CoA (3-(Benzoylamino)-L-Alanine-coenzyme A) is a derivative of coenzyme A. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481