


## 2,3-Dihydroxy-2,3-dihydrobenzoyl-CoA

### Chemical Properties

CAS No. : 809281-20-7  
 Formula: C<sub>28</sub>H<sub>42</sub>N<sub>7</sub>O<sub>19</sub>P<sub>3</sub>S  
 Molecular Weight: 905.66  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



### Biological Description

Description	2,3-Dihydroxy-2,3-dihydrobenzoyl-CoA (2,3-Dihydroxy-2,3-dihydrobenzoyl-coenzyme A) is an aryl-coenzyme A.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1042 mL	5.5208 mL	11.0417 mL
5 mM	0.2208 mL	1.1042 mL	2.2083 mL
10 mM	0.1104 mL	0.5521 mL	1.1042 mL
50 mM	0.0221 mL	0.1104 mL	0.2208 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481