


3-Hydroxy-4-methyloctanoyl-CoA

Chemical Properties

CAS No. : 177261-60-8
 Formula: C₃₀H₅₂N₇O₁₈P₃S
 Molecular Weight: 923.76
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Hydroxy-4-methyloctanoyl-CoA (3-Hydroxy-4-methyloctanoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0825 mL	5.4127 mL	10.8253 mL
5 mM	0.2165 mL	1.0825 mL	2.1651 mL
10 mM	0.1083 mL	0.5413 mL	1.0825 mL
50 mM	0.0217 mL	0.1083 mL	0.2165 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481