

Human Immunoglobulin G4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Human Immunoglobulin G4 (HumanIgG4) is one of the four human IgG subclasses. Unlike other IgG subclasses, IgG4 does not activate complement. In the immune response to filarial parasites, IgG4 may help prevent excessive immune reactions during such infections. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481