

Chitosan-Cy5.5 (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Chitosan-Cy5.5(MW 1000) is a fluorescent dye composed of CY5.5 and chitosan. It is applicable in fields such as bioimaging, drug delivery, and biomaterials research (Ex/Em = 680/710 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481