

Cy5 HA MA (MW 10000)

Chemical Properties

CAS No. :

Formula: 10000 (Average)

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5 HA MA (MW 10000) is a CY5-labeled hyaluronic acid methacrylate widely used in the fields of biomarking and imaging, drug delivery, and tissue engineering and regenerative medicine (Ex/Em = 633/670 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481