


Flagranone C

Chemical Properties

CAS No. :	255064-43-8
Formula:	C ₁₆ H ₁₆ O ₈
Molecular Weight:	336.29
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Flagranone C is a natural product identified within <i>Arthrotrrys flagrans</i> .
Targets(IC50)	Others,Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9736 mL	14.8681 mL	29.7362 mL
5 mM	0.5947 mL	2.9736 mL	5.9472 mL
10 mM	0.2974 mL	1.4868 mL	2.9736 mL
50 mM	0.0595 mL	0.2974 mL	0.5947 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481