


15:0-18:1-15:0 TG

Chemical Properties

CAS No. : 869990-13-6
 Formula: C51H96O6
 Molecular Weight: 805.30
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	15:0-18:1-15:0 TG (1,3-Dipentadecanoyl-2-oleoyl-glycerol) is a type of triglyceride.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2418 mL	6.2089 mL	12.4177 mL
5 mM	0.2484 mL	1.2418 mL	2.4835 mL
10 mM	0.1242 mL	0.6209 mL	1.2418 mL
50 mM	0.0248 mL	0.1242 mL	0.2484 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481