

## 2,3,-Dihydroxybenzoylserine-CoA

### Chemical Properties

CAS No. :

Formula: C<sub>31</sub>H<sub>45</sub>N<sub>8</sub>O<sub>21</sub>P<sub>3</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	2,3-Dihydroxybenzoylserine-CoA (2,3,-Dihydroxybenzoylserine-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481