

AR ligand 41

Chemical Properties

CAS No. :

Formula: C₂₈H₃₅ClN₂O₄

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AR ligand 41 is a PROTAC target protein ligand (Ligands for Target Protein for PROTACs) and can be utilized in the synthesis of BWA-6047.
Targets(IC50)	Ligands for Target Protein for PROTAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481