

## C12 DG PEG (MW 2000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	C12 DG PEG (1,2-Didodecanoyl-rac-glycero-3-methoxypolyethylene glycol) (MW 2000) is a polyethylene glycol lipid that can be utilized in the formation of lipid nanoparticles and liposomes. C12 DG PEG (MW 2000) shows potential for use in cancer research.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481