

(S)-E3 ligase Ligand 42

Chemical Properties

CAS No. :

Formula: C₁₃H₁₀FN₃O₃

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(S)-E3ligaseLigand 42 is an E3 ubiquitin ligase ligand that can connect with a target protein ligand through a linker to form GDC-2992. As a PROTAC, it induces the ubiquitination and degradation of oncogenic proteins.
Targets(IC50)	Ligands for E3 Ligase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481