


Diheteropeptin

Chemical Properties

CAS No. :	202869-98-5
Formula:	C ₂₈ H ₄₂ N ₄ O ₆
Molecular Weight:	530.66
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Diheteropeptin exhibits activity similar to that of transforming growth factor- β and has the ability to inhibit HDAC.
Targets(IC50)	HDAC,Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8844 mL	9.4222 mL	18.8445 mL
5 mM	0.3769 mL	1.8844 mL	3.7689 mL
10 mM	0.1884 mL	0.9422 mL	1.8844 mL
50 mM	0.0377 mL	0.1884 mL	0.3769 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481