

Cy7.5 HA (MW 7000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Cy7.5HA (MW 7000) is a hyaluronic acid labeled with CY7.5, extensively used for cell imaging and in vivo imaging of superficial tissues (Ex/Em = 770/820 nm). |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481