


Abierixin

Chemical Properties

CAS No. :	100634-16-0
Formula:	C40H68O11
Molecular Weight:	724.96
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Abierixin demonstrates significant anticoccidial activity, with mild antibacterial and ionophore properties.
Targets(IC50)	Antibiotic,Parasite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3794 mL	6.8969 mL	13.7939 mL
5 mM	0.2759 mL	1.3794 mL	2.7588 mL
10 mM	0.1379 mL	0.6897 mL	1.3794 mL
50 mM	0.0276 mL	0.1379 mL	0.2759 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481