

U-44069 serinol amide

Chemical Properties

CAS No. :

Formula: C₂₄H₄₁N₅O₅

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	U-44069 serinol amide is a derivative of U-44069 and serves as a stable analogue of the endoperoxide, Prostaglandin H ₂ . It functions as a vasoconstrictor, capable of inducing calcium (Ca ²⁺) influx in preglomerular vessels.
Targets(IC ₅₀)	Prostaglandin Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481