


Fujenal

Chemical Properties

CAS No. :	6750-11-4
Formula:	C ₂₀ H ₂₆ O ₄
Molecular Weight:	330.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fujenal is an abietane-type diterpenoid compound identified in <i>Gibberella fujikuroi</i> . It exhibits moderate inhibitory effects on HeLa tumor cells and serves as an inhibitor of gibberellin biosynthesis. Fujenal holds potential for research in the development of plant growth regulators and anti-tumor agents.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0265 mL	15.1323 mL	30.2645 mL
5 mM	0.6053 mL	3.0265 mL	6.0529 mL
10 mM	0.3026 mL	1.5132 mL	3.0265 mL
50 mM	0.0605 mL	0.3026 mL	0.6053 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481