


(R)-3-Hydroxy-5-phenylvaleryl-CoA

Chemical Properties

CAS No. : 473994-63-7
 Formula: C₃₂H₄₈N₇O₁₈P₃S
 Molecular Weight: 943.75
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|---|
| Description | (R)-3-Hydroxy-5-phenylvaleryl-CoA ((R)-3-Hydroxy-5-phenylvaleryl coenzyme A) is a derivative of coenzyme A. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 1.0596 mL | 5.298 mL | 10.596 mL |
| 5 mM | 0.2119 mL | 1.0596 mL | 2.1192 mL |
| 10 mM | 0.106 mL | 0.5298 mL | 1.0596 mL |
| 50 mM | 0.0212 mL | 0.106 mL | 0.2119 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481