


L-2,5-Dihydrophenylalanine

Chemical Properties

CAS No. :	16055-12-2
Formula:	C9H13NO2
Molecular Weight:	167.21
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	L-2,5-Dihydrophenylalanine can induce apoptosis in human medulloblastoma leukemia cells [HL-60].
Targets(IC50)	Apoptosis, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9805 mL	29.9025 mL	59.805 mL
5 mM	1.1961 mL	5.9805 mL	11.961 mL
10 mM	0.5981 mL	2.9903 mL	5.9805 mL
50 mM	0.1196 mL	0.5981 mL	1.1961 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481