

BMS-986497

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BMS-986497 (ORM-6151) is a CD33-targeted antibody conjugated with a GSPT1 degrader. It delivers the GSPT1 degrader SMol006 to CD33-expressing cells, inducing the degradation of the GSPT1 protein. BMS-986497 holds potential for research in acute myeloid leukemia (AML).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481