

LSN3172176 Precursor

Chemical Properties


CAS No. :

Formula: C₂₁H₂₈BrN₃O₃

Molecular Weight: 450.38

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	LSN3172176 Precursor is a precursor for ¹¹ C-radiolabelling LSN3172176 PET scan.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2203 mL	11.1017 mL	22.2035 mL
5 mM	0.4441 mL	2.2203 mL	4.4407 mL
10 mM	0.222 mL	1.1102 mL	2.2203 mL
50 mM	0.0444 mL	0.222 mL	0.4441 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481