

(2-Pyridyldithio)-PEG4-alcohol

Chemical Properties

CAS No. : 851961-99-4


Formula: C13H21NO4S2

Molecular Weight: 319.44

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2-Pyridyldithio)-PEG4-alcohol, a PEG-based PROTAC linker, is used for the synthesis of PROTACs[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker: one binds to an E3 ubiquitin ligase, and the other targets a specific protein. They harness the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1305 mL	15.6524 mL	31.3048 mL
5 mM	0.6261 mL	3.1305 mL	6.261 mL
10 mM	0.313 mL	1.5652 mL	3.1305 mL
50 mM	0.0626 mL	0.313 mL	0.6261 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481