

1× PBS Buffer (Sterile)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	1× PBS Buffer (Sterile), namely phosphate-buffered saline, is a commonly used buffered salt solution exhibiting excellent pH buffering capacity and isotonic properties. It is suitable for cell culture, biochemical experiments, and molecular biology research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481