


## Tetrahydrocorticosterone

## Chemical Properties

CAS No. :	68-42-8
Formula:	C <sub>21</sub> H <sub>34</sub> O <sub>4</sub>
Molecular Weight:	350.49
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tetrahydrocorticosterone (THB) is a glucocorticoid (GC) that has anti-inflammatory functions. For instance, 5 $\alpha$ THB inhibits inflammatory actions by binding to GC receptor. HPLC studies in rat brains have reported that tetrahydrocorticosterone (THB) isomers are differentially distributed in males and females.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8531 mL	14.2657 mL	28.5315 mL
5 mM	0.5706 mL	2.8531 mL	5.7063 mL
10 mM	0.2853 mL	1.4266 mL	2.8531 mL
50 mM	0.0571 mL	0.2853 mL	0.5706 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481