

Ascorbic Acid Assay Kit (UV Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Ascorbic Acid (AsA), among other functions, are also known as Vitamin C, acts as a coenzyme, free radical scavenger, electron donor/acceptor, and substrate for the biosynthesis of oxalate and tartrate. As the most important antioxidant in plant cells, AsA plays a crucial role in protecting chloroplasts from oxidative damage and is also an important indicator for measuring the quality of agricultural products.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481