

## Asparaginase Assay Kit (Spectrophotometry)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Asparagine synthase (AS) is a type of aminotransferase widely found in living organisms, which catalyzes the transfer of an amino group from glutamine to aspartic acid; the formation of asparagine is a detoxification reaction when plants are under ammonia toxicity.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481