

Blood Magnesium Concentration Assay Kit (Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Magnesium is an activator for various enzymes, such as phosphatases, creatine kinase, hexokinase, and carboxylases. Magnesium is also an essential element for the structure of macromolecules like DNA, RNA, and ribosomes. It is a crucial element for maintaining normal nerve and muscle function. Deviations of serum magnesium concentration from the normal range are associated with certain kidney and endocrine diseases.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481