

CAD Assay Kit (Microanalysis)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Cinnamyl-alcohol dehydrogenase (CAD) is one of the key enzymes in the lignin biosynthetic pathway, catalyzing the final step in the synthesis of lignin monomers. It catalyzes the reduction of various cinnamaldehydes (such as coumaraldehyde, sinapaldehyde, and coniferaldehyde) to their corresponding cinnamyl alcohols. This enzyme is predominantly found in higher plants, yeast, and fungi. Research on this enzyme can help elucidate the metabolic mechanisms of lignin deposition during the development of various biological cells, providing a basis for reducing stone cell content in fruits and improving their quality.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481