

CHI Assay Kit (Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Chalcone isomerase (CHI) is the first recognized enzyme involved in the synthesis of flavonoids and is one of the key enzymes in the flavonoid metabolic pathway. Together with chalcone synthase, it constitutes the rate-limiting enzyme for the biosynthesis of flavonoids.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481